

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
	AD-A123139	
4. TITLE (and Subtitle) A Table for Rapid Determination of Log ₁₀ LD ₅₀ or ID ₅₀ Values		5. TYPE OF REPORT & PERIOD COVERED N/A
		6. PERFORMING ORG. REPORT NUMBER
7. AUTHOR(s) Ralph W. Kuehne		8. CONTRACT OR GRANT NUMBER(s)
9. PERFORMING ORGANIZATION NAME AND ADDRESS U.S. Army Medical Research Institute of Infectious Diseases SGRD-UIZ-S Fort Detrick, Frederick, MD 21701		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS N/A
11. CONTROLLING OFFICE NAME AND ADDRESS U.S. Army Medical Research and Development Command, Office of the Surgeon General Dept. of the Army, Washington, DC 20314		12. REPORT DATE 7 Oct 82
		13. NUMBER OF PAGES 6 pages
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		15. SECURITY CLASS. (of this report) UNCLASSIFIED
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for public release; distribution unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES To be published in the Journal of Clinical Microbiology		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) LD ₅₀ , 50% endpoint		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) A reference table is presented for rapid determination of LD ₅₀ or ID ₅₀ values when administering 10-fold serial dilutions to six animals per dilution.		

ADA123139

DTIC
ELECTE
JAN 10 1983

S D

DTIC 123139

FORM 1073 JAN 73

EDITION OF 1 NOV 65 IS OBSOLETE

88 01 10

UNCLASSIFIED
SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

Accession For	
NTIS GRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability Codes	
Dist	Avail and/or Special
A	



A Table for Rapid Determination of Log_{10} LD_{50} or ID_{50} Values

RALPH W. KUEHNE

U.S. Army Medical Research Institute of Infectious Diseases, Fort Detrick, Frederick,
Maryland 21701

The views of the author do not purport to reflect the positions of the Department
of the Army or the Department of Defense

Phone No. (301) 663-7373

"Approved for public release; distribution unlimited"
Cleared on 30 Sept 82

83 01 10 071

ABSTRACT

A reference table is presented for rapid determination of LD₅₀⁷ or ID₅₀^M values when administering 10-fold serial dilutions to six animals per dilution.

While the determination of an LD₅₀ or ID₅₀ is neither a complicated nor a difficult procedure (1), many laboratories continue to perform these calculations on each titration manually. In some of our laboratory areas where access to a computer is not readily available, the table (Table 1) presented below has been of great value in saving time, reducing errors, and enabling even untrained personnel to record LD₅₀ values from titration results. Values shown are for titrations utilizing six animals per point, a commonly used number for titrating viruses or sera for antibody (2). The method of Reed and Muench (1) "utilizes, at least in part, the large, total number of animals involved. The effect is that of using, at the two critical dilutions between which the endpoint lies, larger groups of animals than were actually included at these dilutions."

To use the table, let A represent the log₁₀ of the highest dilution showing the desired effect (mortality, infectivity, immunity, survival, etc.) in 6/6 animals, and B through E the log₁₀ of the four next consecutive 10-fold dilutions.

Typical titration examples are given in Table 2. Example number 2 shows a titration giving the following results: $10^{-2} = 6/6$, $10^{-3} = 6/6$, $10^{-4} = 4/6$, $10^{-5} = 2/6$, $10^{-6} = 2/6$ and $10^{-7} = 0/6$. In this instance, A = 3 (the log₁₀ of the highest dilution showing the effect in all animals), B = 4, C = 5, etc. The value from Table 1 would be derived from the line showing 6, 4, 2, 2, 0 (bracketed results in Table 2), or B.75. Since B = 4 in this example, the log₁₀ LD₅₀ would equal 4.75.

The table assumes that no more than five dilutions are required to accomplish 100 to 0% response, and that an orderly dose-response relationship exists. Since the mantissa of each value shown represents the proportionate distance between two adjacent dilutions, this number can also be applied to calculate the 50% endpoint when dilutions other than 10-fold are used, i.e., to other log bases.

LITERATURE CITED

1. Reed, L. J., and H. Muench. 1938. A simple method of estimating fifty per cent endpoints. Am. J. Hyg. 27:493-497.
2. Lennette, E.H. 1964. General principles underlying laboratory diagnosis of viral and rickettsial infections, p. 1-66. In E. H. Lennette and N. J. Schmidt (ed.), Diagnostic procedures for viral and rickettsial diseases, 3rd ed. American Public Health Association, New York.

TABLE 1. Log_{10} LD₅₀ or ID₅₀ Values for 6 Animals/Dilution

*10 ^{-A} 10 ^{-B} 10 ^{-C} 10 ^{-D} 10 ^{-E} LD ₅₀						*10 ^{-A} 10 ^{-B} 10 ^{-C} 10 ^{-D} 10 ^{-E} LD ₅₀						
6	5	5	5	5	0	D.20	6	4	3	1	0	B.84
6	5	5	4	0	0	D.00	6	4	3	0	0	B.69
6	5	5	3	0	0	C.71	6	4	2	2	0	B.75
6	5	5	2	0	0	C.53	6	4	2	1	0	B.63
6	5	5	1	0	0	C.40	6	4	2	0	0	B.50
6	5	5	0	0	0	C.30	6	4	1	1	0	B.47
6	5	4	4	0	0	C.80	6	4	1	0	0	B.36
6	5	4	3	0	0	C.55	6	4	0	0	0	B.25
6	5	4	2	0	0	C.38	6	3	3	3	0	C.00
6	5	4	1	0	0	C.24	6	3	3	2	0	B.83
6	5	4	0	0	0	C.13	6	3	3	1	0	B.67
6	5	3	3	0	0	C.33	6	3	3	0	0	B.50
6	5	3	2	0	0	C.16	6	3	2	2	0	B.59
6	5	3	1	0	0	C.00	6	3	2	1	0	B.45
6	5	3	0	0	0	B.84	6	3	2	0	0	B.31
6	5	2	2	0	0	B.88	6	3	1	1	0	B.29
6	5	2	1	0	0	B.76	6	3	1	0	0	B.16
6	5	2	0	0	0	B.64	6	3	0	0	0	B.00
6	5	1	1	0	0	B.60	6	2	2	2	0	B.38
6	5	1	0	0	0	B.50	6	2	2	1	0	B.20
6	5	0	0	0	0	B.40	6	2	2	0	0	B.00
6	4	4	4	0	0	C.63	6	2	1	1	0	B.00
6	4	4	3	0	0	C.41	6	2	1	0	0	A.88
6	4	4	2	0	0	C.25	6	2	0	0	0	A.75
6	4	4	1	0	0	C.12	6	1	1	1	0	A.80
6	4	4	0	0	0	C.00	6	1	1	0	0	A.70
6	4	3	3	0	0	C.17	6	1	0	0	0	A.60
6	4	3	2	0	0	C.00	6	0	0	0	0	A.50

* Serial Log Dilutions.

TABLE 2. Examples of Typical Titration Results

Example	10^{-1}	10^{-2}	10^{-3}	10^{-4}	10^{-5}	10^{-6}	10^{-7}	10^{-8}	Value of A	LD ₅₀ Value From Table 1	Log ₁₀ LD ₅₀
1.	6/6	6/6 ^A	5/6 ^A	4/6 ^C	2/6 ^D	0/6 ^E	0/6		2	C.38	4.38
2.		6/6	6/6 ^A	4/6 ^B	2/6 ^C	2/6 ^D	0/6 ^E		3	B.75	4.75
3.		6/6	6/6	5/6 ^A	4/6 ^B	3/6 ^C	1/6 ^D	0/6 ^E	4	B.84	5.84

* 10-Fold Dilutions.